BLOK-JET

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BLOK-JET

COMPRESSED AIR AND DENTAL SURGERY ASPIRATION PRODUCTION PLANTS

1) BESPOKE DESIGN

This first section (to paragraph 4 included) highlights Blok-Jet's most innovative features and is followed by instructions on its installation, functioning and servicing.

Blok-Jet is a system that allows you to put together a compact and centralized compressed air and dental suction plant of any size tailored to your requirements.

Installing several heads and suction units reduces the possibility of a total stoppage in the compressed air and suction supply and also saves energy in the motor start-up phase.

2) PROFESSIONAL PERFORMANCE

The introduction of the inverter along with the software programme has allowed us to create systems offering a level of professional performance that was until now considered unachievable.

The inverter (also known technically as the VSD*) is an electronic drive that optimizes the functioning of electric motors; the software is an artificial intelligence that constantly monitors the functioning of the machines as well as the plant's overall operation.

The VSD and the software programme reduce stress on the motors and automatically keep air-flow and vacuum levels constant and adjust the capacity to meet demand in real time.


* VSD is included only in the suction units and the software and remote control are present both in the suction units and in the compressed air units.

3) SAFETY AND PREVENTION

Where an anomaly arises such as overload, high temperature, power fluctuation or the risk of flooding, **the system goes into self-protection mode without stopping and adjusts the plant's operation to emergency mode.** Once the anomaly has passed, the plant automatically returns to normal operation.

Also, if the system's safety is jeopardized, **the programme sets off an alarm** that can be transferred to a supervised location or, **via the internet**, to the technician on duty; at the same time, the plant's display screen shows an image of the faulty machine so that the technician on duty can decide whether to **intervene remotely.**

4) ENERGY SAVING AND ENVIRONMENTAL SUSTAINABILITY

The suction curves shown highlight the energy saving made when using a Uni-Jet 501 equipped with VSD compared to a Uni-Jet 1000 equipped with a supplementary valve: **at a capacity of 216 m³/h, the energy saving is 7.5 kW /h; at 480 m³/h the saving is 3 kW/h.** At 585 m³/h, the power absorption of the two suction units is almost the same.

As well as saving energy, there is also a saving in raw materials: **for a clinic composed of 32 dental units, a 501 suction group weighing 88 kg is sufficient while for the same number of dental units, without VSD, a Uni-Jet 1000 weighing 155 kg is required.** The first has a power capacity of 7.5 kW and the second a power capacity of 15 kW.


This is not just a centralized plant but one that has been planned, designed and manufactured to meet professional standards as well as to achieve significant savings in energy and raw materials.

5) WHY USE AN OIL-FREE COMPRESSOR

The old system of oil-bath rotary compressors is still used in large installations as it was believed these were needed to supply the large air-flow demand required by the number of positions in hospital installations. Now we know that **dentistry requires clean, dry and filtered air, both for delicate professional operations such as surgery and prosthetic work as well as for maintaining the pneumatic equipment used in dentistry in good working order.** In lubricated compressors, the oil heats up when it is in operation and, despite the filters, the finest of particles still pass in the air through all possible barriers compromising fillings and risking the cementation of the prosthesis.

The need for dry, filtered and sterile air made us determined to solve this problem and our centralized air compression units are safe, replacing rotary oil-bath compressors of any capacity with the added benefit of notable savings in energy and maintenance costs.

Cattani S.p.A. can demonstrate how an oil-free compressor can have the same, or longer, longevity than lubricated compressors: it is a question of choosing the right materials and of manufacturing precision.

1a) PLANT FUNCTIONING

When the compressed air centralized unit is switched on, the heads will start one by one in sequence a few seconds apart and will remain running until the maximum set pressure is reached; they will re-start in the same way when the pressure drops to the set minimum level. The heads draw in air from the environment where the compressor is installed. On request, a connection is available to supply the heads with air from outside the plant room in order to draw in cleaner and fresher air. By channelling the air drawn in, the sound vibration can be reduced by one or two dB (A) on one- or two-cylinder heads and three dB(A) on three-cylinder heads.


Each head is equipped with a one-way valve; a head stopping does not interrupt the Blok-Jet's functioning as the one-way valve is always closed when the head is stopped and the plant will continue to function normally with one less head. The head stopping is shown on the control box display; the same message can be transferred to another location via the internet.

2b) COMPRESSED AIR TREATMENT

During compression, air heats to a temperature of up to 100°C; on exiting the heads, the air is transferred to the air/air exchanger where it is brought down to ambient temperature + 5 °C. The air then passes to the centrifugal separator where the water separated during the drop in temperature is deposited. After the first drying process, the air passes through the pre-filter where it is cleaned of solid particles and dust (99.99% filtration at 3 micron). The filtered air is then transferred to the refrigerated dryer where it deposits any residual humidity (ISO 4 dryness classification - dew point rating of 3°C under pressure). From this last dryer, the air passes through the post-filter (filtration 99.99999% at 0.1 micron) and is then transferred to the tank ready for use.

As an alternative to the refrigerated cycle described above, an air drying system is available which reaches a dew point adjustable between -20 °C a -40°C.

STANDARD PLANT

3c) FUNCTIONING OF SUCTION UNITS

When there are few operators working, only one suction unit will operate at low revolutions and all the suction piping will quickly reach the set level. If the demand increases, the software will activate the inverter to increase the Hz available and therefore the rotation speed of the first machine; once the set maximum speed of the first aspirator is reached, the second suction unit will be activated and both machines will function at a reduced speed. The system will continue to adjust the rotation speed and activate or deactivate the available machines according to constant suction level set by the programme and using the least possible amount of energy.

4d) FUNCTIONING OF THE SEPARATOR TANK

The aspirated fluid transferred from operating theatres (green route) is filtered before entering the separator tank. Inside the tank, the air (blue route) separates from the liquids (red route); the liquid accumulates at the bottom of the tank while the air exits at the top. When the liquid reaches the maximum set level, the level control (driven by a set of three probes) activates the Turbo-Smart as draining pump; the pump stops when the minimum set level is reached. The entire operation is repeated continuously and automatically when the plant is in operation without inhibiting suction. Air exiting from the separator tank joins the air exiting from the Turbo-Smart and follows the route to the anti-bacterial filter where it is transferred to the two suction groups and expelled outdoors.

5e) ISO 18 HYDROCYCLONE (ISO 11143 AMALGAM SEPARATOR)

The liquid coming from the Turbo-Smart (as draining pump) rotates down the hydrocyclone and reverses direction and rotation when it reaches the narrowest section of the cone: at this point the particles of amalgam are deposited on the bottom of the tank connected to the hydrocyclone's cone while the liquid, cleaned of amalgam, exits at the top of the cone. The ISO hydrocyclone has a capacity of 18/l min.



The system does not start to drain immediately; when the liquid has reached separation speed the drainage valve will open automatically. The ISO 18 hydrocyclone separates 98.5% of the amalgam sample contained in the test measure established by the ISO standard.

Both an audible and visible alarm will be activated when the level of amalgam in the collection tank has reached 90%. The audible alarm may be switched off and the system restarted but the light will remain on. After around one week's work, the same alarms will signal that the tank has reached 100% capacity and the system will stop; the amalgam collection tank must be emptied before the system can be restarted. Before closing the container, it is recommended that the empty space be filled with an authorized disinfectant, then sealed and given to authorised disposal experts as regulated by the relevant legislation in force.

THE PLANT ROOM

The plant room should be prepared before the machines arrive: installing air compression and suction units in a construction site is often difficult and not without unexpected factors. It is recommended to provide for an external air feed for the compressor as well as a floor drain and a vent for the air exhausted from the suction groups, preferably to the roof.

Machines create heat when operational and an environmental temperature higher than 40°C may damage the motors and compromise the drying process for the compressed air.

Plant rooms without windows are extreme but not rare; however, windows are often not enough to maintain an acceptable room temperature. It is therefore advised to install an air exchanger or air conditioning in the plant room to ensure a minimum room temperature of + 5 °C and a maximum of + 35°C. The conditioned air must flow directly onto the machines. It is also advised to use the expertise of an air-conditioning technician for the planning and installation of any temperature control equipment.

The power cables, control boxes and machines must be protected against direct or indirect electrical contact, and from power overloads or over-current pursuant to Regulations CEI 64-8 for class one equipment. Correspondent IEC data is available on request.

INSTALLATION

It is recommended to the person(s) installing the plant to consult the instruction manual and to contact the manufacturer for any further clarification or requirement.

All pipes connecting the machines together and to other installations must always be separated by a small section of flexible tubing to dampen vibrations caused by moving parts. The flexible connections should be positioned only in the plant room and where they can be easily inspected.

The various sections of rigid piping for the suction canalization should be welded together using only glues and systems recommended by the manufacturer of the canalization used.

The compressed air plant should be connected using a reinforced flexible hose compliant with EN 853 Standard.

Do not twist any pipes and check that they are all air-tight: any leaks in compressed air or suction (blue route), will inhibit the correct functioning of the machines as well as the compressed air and suction network.

Air and liquid from operating theatres must always be treated as contaminated.

Check that the tubes are clean before switching on any machines; heavy debris, as well as water or dust, may damage the machines and harm operators. Check the sense of rotation of the suction units, compressor heads, fans and drainage pumps.

TESTING

Not all dentists' units are the same and the suction requirements may vary according to the diameter of external and internal tubes, the number of bends or size of the selection valves. There may be different results based on the models of dentist's units present in the same installation. The actual testing can only therefore be carried out on the centralized production plants.

Procedure: Compressor charging times and suction delivery head should be tested at all chairs if the installation was intended for the simultaneous use of all chairs; the centralized unit should be able to supply 45/50 l/min. of compressed air to each chair at a minimum pressure of 3 bar, while suction should be maintained at the set level of 250/300 litres to each chair.

The working pressure of compressed air and the delivery head should then be adjusted as required by the operators. Please note that the pressure and delivery head levels have a tolerance due to the automatic intervention of the motors; this tolerance does not inconvenience operators during their operations.

It is advised to check the correct functioning of the compressed air and suction units, in particular: one-way valves, supplementary valve, exhausted air vent, liquids drainage, separator tank levels, motor fans, air/air exchanger ventilator on the compressor and the silica gel regeneration.

Particular attention should be paid to checking the programming of the control panels, from the operating functions to worker safety, as well as the operating or danger signalling and the automatic safety devices on the machines themselves.

MAINTENANCE OF THE COMPRESSED AIR AND SUCTION CENTRALIZED PLANT

We advise those responsible in the clinics or surgeries to engage a qualified technician to carry out the servicing and to use only original spare parts. Only six-monthly checks are required on compressors and suction units that are located in an environment that is clean and where the temperature is controlled (from + 5° C to + 35° C). The technician should carry out all the checks listed for the inspection and which are repeat below.

Check the compressor charging times and suction delivery heads at all chairs. Check the correct functioning of the compressed air and suction units, in particular: one-way valves, supplementary valve, exhausted air vent, liquids drainage and separator tank levels. Check the correct functioning of the motor fans, air/air exchanger fan on the compressor and the condition of the silica gel regeneration.

Particular attention should be paid to check the programming of the control panels, the operating functions, worker safety, operational or danger signalling, and automatic safety interventions on the machines themselves.

Pipes, seals and membranes should be replaced each time they are disassembled for interventions of any nature, whether located in a temperature-controlled environment or not; it is recommended that they be replaced every 12/18 months. Motor bearings should be replaced every ten thousand hours or every five working years. Anti-bacterial and absolute filters should be replaced once a year; filters for aspirated secretions should be checked monthly, more often if the filter on the dentist's chair has a wider gauge mesh.

The technician will use his or her experience as well and look at the amount of work under which the units are placed, the noise made by the machines, their temperature when in operation and any air or liquid leakages to decide whether to shorten the intervals between services and whether or not to replace any parts.

Where machines are installed outdoors, check that there are no obstacles to the free ventilation of the machines (paper, leaves or other).

The maintenance technician must not alter the machines or their functioning, or the safety devices in any way; in particular, no welding must be carried out on the compressor tank.

CLEANING AND DISINFECTION OF THE DENTAL SURGERY ASPIRATION EQUIPMENT

The following paragraph is aimed mainly at surgery staff in charge of disinfection of the equipment. At night, at the end of the working day (more often if necessary), suck a detergent and disinfectant solution (Puli-Jet plus with anti-scale agent) directly on the dentist's units using the Pulse Cleaner; the pulsation in fact creates the turbulence necessary to dislodge dirt so that the disinfectant solution can work properly. Continuous use of Puli-Jet Plus New with anti-scale agent dissolves blood and mucus and dislodges scale from the aspirator walls and the canalization thereby removing the anchor for the bacterial biofilm. Puli-Jet Plus New with anti-scale agent can remain inside the aspirator and the canalization as long as required without damaging the surfaces being treated. Therefore, night-time disinfection is advised as well as before holidays or weekends; it is also recommended not to rinse the equipment after disinfection.

Often, a dentist's surgery works late into the night then there is no more time. But, disinfection carried out in the morning, a few minutes before starting work, does not have the same affect as overnight disinfection: treatment the day after cannot cancel the biological changes that occurred during the night in the liquids suctioned the day before. Blood and mucus, agitated by forced air, generate a large quantity of foam; after treatment with disinfectant detergents and daily cleaning of the secretions filter on the dentist's unit, it is **advised to also place a anti-foam disinfectant tablet in the actual filter.**

The day after, **the liquid aspirated through the filter will pick up some particles of anti-foam and disinfectant which will inhibit the formation of foam as well as deliver a continuous anti-bacterial action;** this is better than using more aggressive products that may damage machinery and equipment.

Methodical cleaning of the aspirator will ensure its proper functioning while a lack of cleaning or disinfection will promote the formation of layers of bacteria. Over time, these will reduce suction leading to moving parts overheating and seizing up.

CONCLUSIONS

We have endeavoured to use modern technology to build professional centralized plants that are safe, flexible and adaptable offering the possibility of remote control and adjustment through the internet or wireless connection.

During our research and the manufacture of the units, we realised that using modern technology inevitably led us to save energy and raw materials, savings that mankind needs to make if it wants to protect our planet.

We warmly thank all those* who have already shown their appreciation for our commitment and those who will do so in the future.

Cattani S.p.A.

* There are really many - professionals, constructors, retailers and technicians.

CATTANI S.P.A. IN THE WORLD

AUTHORISED DEALERS ABROAD

CATTANI AUSTRALIA	andrew@cattani.com.au bill@cattani.com.au	Australia and New Zealand
DENT UNIT	kadlecek@dentunit.cz	Czech Republic
DEN-TEC APS	info@den-tec.dk	Denmark
FONTAIN HERVE' S.A.R.L.	herve@fontoin.fr	France
CATTANI DEUTSCHLAND GMBH & CO. KG	info@cattani.de	Germany
DENT-EAST CO.LTD.	Iren.Fuzesi@dent-east.com	Hungary
HOODIANN CO.	mercede.naderi@hoodian.com	Iran
DENTAL SOR A.S.	post@dentalnet.no	Norway
CATTANI POLSKA	office@cattani-polska.com.pl	Poland
SITEA ROMANIA SRL	fabrizio@sitea.ro	Romania
STOMATORG	e.devletbaeva@stomatorg.ru	Russia
CORAL LTD	international.coral@gmail.com	Russia
AL-AJAJI TRADING COMPANY	eng.mahmoud@ajajigroup.com	Saudi Arabia - KSA
HEUNG JIN MEDICAL CO.LTD	heungjinmd@yahoo.co.kr	South Korea
INCOTRADING S.A.	stig@incotrading.net andres@incotrading.net	Spain, Portugal, Central and South America
GUNEY DYS DEPOSU SAN. TIC.A.S	Cansun.Mert@guneydis.com	Turkey
AL-HAYAT PHARMACEUTICALS	m.abueida@alhayatdental.com	United Arab Emirates - UAE
CATTANI ESAM U.K. LIMITED	Steve@cattaniesam.co.uk	United Kingdom